

## **Human molecular genetics: information for sampling and where to send**

### **Molecular studies:**

Genomic DNA is extracted from EDTA whole blood (3 ml) and the concentration is adjusted to a range of 25–50 ng/μl.

All samples should be labeled with a sticker decal with patient's name and date of birth.

All samples must be sent with a requisition form that includes the clinical history and functional parameters essential for the interpretation of the laboratory results (Annex 3).

### **Analysis request:**

This is available in the Haemostasis and Molecular Haematology Laboratory of Centro Hospitalar e Universitário de Coimbra. All information (Specimen Collection, Handling and Transportation) can be found at [www.chuc-hematologiaclinica.org](http://www.chuc-hematologiaclinica.org)